

<http://lifescience.fujifilm.com>

QuickGene Video  
"All about QuickGene"



**NOW SHOWING!**

Application Guide  
"No.2-49"



**Available on our website**

FUJIFILM

Life Science

# Nucleic Acid Isolation System Selection Guide

For QuickGene Selection Guide



FUJIFILM

FUJIFILM Corporation 7-3, Akasaka 9-Chome, Minato-ku, Tokyo 107-0052, Japan, Tel:+81-3-6271-2158, Fax:+81-3-6271-3136, E-mail: sginfo@fujifilm.co.jp

FUJIFILM Europe GmbH Heesenstr. 31, 40549 Dusseldorf, Germany, Tel: +49-211-5089-144, Fax: +49-211-5089-9144, E-mail: lifescience@fujifilm.eu

FUJIFILM UK Ltd. Unit 12 St Martins Way, St Martins Business Centre, Bedford, MK42 0LF, U.K, Tel: +44-1234-572226, Fax: +44-1234-245293, E-mail: lifesciences@fuji.co.uk

富士胶片(中国)投资有限公司 31st floor, Hong Kong New World Tower, No.300 Huai Hai Zhong Road, Shanghai, P.R China, Tel: +86-21-3302-4655 ext.363, Fax: +86-21-6384-3322, E-mail: wgxian@fujifilm.com.cn

[http:// lifescience.fujifilm.com](http://lifescience.fujifilm.com)

Specifications and system configuration subject to change for improvement without notice. All other product names mentioned herein are the trademarks of their respective owners.

宣-08-5-(F1099)・3-3(F1099) 宣伝物コード6007 6786

2008.5

# QuickGene Series

## Covers a wide range of areas to realize your ideas.

The "QuickGene" series uses a porous membrane developed through the application of FUJIFILM's membrane production technology to realize high purity and high yield in nucleic acid extraction. Versatile extraction kits support various samples to expand the application and possibility of DNA / RNA extraction, from basic research to medicine, food, agriculture and forensic criminal investigations.

### Extraction kits features

### Quick and easy DNA / RNA extraction with QuickGene kits

#### All-In-One Package

Sample preparation can be conducted with the reagents, enzymes and vessels included in a single package. Nucleic acid extraction can be conducted as soon as the kits arrive.

#### Store at room temperature

Store the reagents at 15°C ~ 28°C. No need for refrigerated storage.  
\*For enzyme reagents, refrigerated storage is recommended after use.

#### No hazardous organic solvents


The cartridges and solvents are all supplied without DNase and RNase to avoid contamination. Environmentally friendly extraction can be conducted without using hazardous organic solvents.

#### Compact size

To minimize space requirements, all necessary items are packaged in a single compact package. Kit S for QuickGene-810 / Mini80, SP kit contains 96 samples, and kit L for QuickGene-610 contains 48 samples.

### Core technology for high-purity and high-yield extraction

The nucleic acid adsorptive medium used in QuickGene series is a porous membrane developed through application of FUJIFILM's advanced polymer membrane production technology. It is only 80 μm thick, making it incomparably thinner than conventional glass fibers. Because of the outstanding adsorptive and desorptive performances of the membrane, nucleic acid can be rapidly and reliably extracted at low pressure without being damaged, which realizes high-quality nucleic acid extraction.



## One for each person QuickGene-Mini80



### Features

The series' smallest system enabling nucleic acid extraction through simple operation; just set the sample and rotate the pressurizing switch. No need to move from the lab bench throughout the extraction.

### Features

compact reasonably priced

#### Extraction kits (seven)

DNA	whole blood; tissue; plasmid II
RNA	blood cell; tissue II; cultured cell; cultured cell HC

#### Specifications

##### Overview

- Throughput: 1 to 8 samples per run

##### Operating conditions

- Supply voltage: AC 100-240 V
- Power supply frequency: 50/60 Hz
- Temperature: 15-30°C
- Humidity: 30-80 % (non-condensing)

##### Physical specifications

- Dimensions: 280 (W) x 220 (D) x 180 (H) mm
- Weight: Approx. 3 kg

## Desktop multifunction model QuickGene-810



### Features

A multifunctional automated system realizing high-purity high-yield DNA / RNA extraction from varied samples (human, mouse, wheat, E.coli, cell, etc.).

### Features

automated multifunctional

#### Extraction kits (seven)

DNA	whole blood; tissue; plasmid II
RNA	blood cell; tissue II; cultured cell; cultured cell HC

#### Specifications

##### Overview

- Automated stages: Sample binding, washing and elution
- Throughput: 1 to 8 samples per run
- Display: LCD (16 characters x 1 line)

##### Operating conditions

- Supply voltage: AC 100-240 V
- Power supply frequency: 50/60 Hz
- Temperature: 15-30°C
- Humidity: 30-80 % (non-condensing)
- Power Consumption: 65 W

##### Physical specifications

- Dimensions: 450 (W) x 330 (D) x 400 (H) mm
- Weight: Approx. 21 kg

##### Options

- Carriage sets
- Sample tube rack (sold in sets of 4)

## Series' largest extraction scale QuickGene-610L



### Features

A stable high-purity high-yield extraction system, enabling automated extraction of approx. 50μg DNA from 2 ml whole blood sample\*. Suitable for checking multiple parameters using limited amounts of blood in clinical research or livestock /animal research.

### Features

automated large-scale

#### Extraction kits (one)

DNA	whole blood L
-----	---------------

#### Specifications

##### Overview

- Automated stages: Sample binding, washing and elution
- Throughput: 1 to 6 samples per run
- Display: LCD (16 characters x 1 line)

##### Operating conditions

- Supply voltage: AC 100-240 V
- Power supply frequency: 50/60 Hz
- Temperature: 15-30°C
- Humidity: 30-80 % (non-condensing)
- Power Consumption: 100 W

##### Physical specifications

- Dimensions: 580 (W) x 330 (D) x 400 (H) mm
- Weight: Approx. 24 kg

## Spin-cartridge multifunctional kit QuickGene SP kit



### Features

Rapid and easy DNA / RNA extraction using equipment already available in your laboratory, such as centrifuges and microtubes. Because washing and recovery of nucleic acid can be performed in a tabletop compact centrifuge, work efficiency can be dramatically improved.

### Features

Spin method

#### Extraction kits (six)

DNA	whole blood (spin method); tissue (spin method); plasmid II (spin method)
RNA	tissue (spin method); cultured cell (spin method); cultured cell HC (spin method)

# DNA Kit

For QuickGene-810/Mini80 [for 96 samples]

## DNA tissue kit

Processing time: 13 min/ 8 samples

Pretreatment enzyme, Tissue lysis buffer, Lysis buffer, Wash buffer, Elution buffer, Cartridges, Caps, Collection tubes, Waste tubes

## Plasmid kit II

Processing time: 6 min/ 8 samples

Pretreatment enzyme, Lysis buffer, Resuspension buffer, Alkaline solution, Neutralization buffer, Wash buffer, Elution buffer, Cartridges, Caps, Collection tubes, Waste tubes

## DNA whole blood kit

Processing time: 6 min/ 8 samples

Pretreatment enzyme, Lysis buffer, Wash buffer, Elution buffer, Cartridges, Caps, Collection tubes, Waste tubes

For QuickGene-610L [for 48 samples]

## DNA whole blood kit L

Processing time: 12 min/ 6 samples

Pretreatment enzyme, Lysis buffer, Wash buffer, Elution buffer, Cartridges, Waste tubes

For spin-cartridge method extraction [for 96 samples]

## DNA tissue kit (spin method)

Pretreatment enzyme, Tissue lysis buffer, Lysis buffer, Wash buffer, Elution buffer, Cartridges, Waste tubes

## Plasmid kit II (spin method)

Pretreatment enzyme, Lysis buffer, Resuspension buffer, Alkaline solution, Neutralization buffer, Wash buffer, Elution buffer, Cartridges, Waste tubes

## DNA whole blood kit (spin method)

Pretreatment enzyme, Lysis buffer, Wash buffer, Elution buffer, Cartridges, Waste tubes

Extraction example

ca. 4 µg/ 5 mg Balb/c Mouse tail

ca. 12.5 µg/ 1 ml culture/GAPDH/ DH5α

ca. 5 µg/ Whole blood 200 µl

ca. 50 µg/ Whole blood 2 ml

### Mammalian (Human/Cow/Poultry/Dog/Cat)

- DNA isolation for genetic test
- Genotyping
- Identification and genotyping

Nail	16
Dental pulp and hard tissues (teeth and bones)	17
Paraffin-embedded samples	46
Oral swab	-

Hair	-
Lymphatic node, Liver, Kidney	-
Blood Spot	-

Whole blood	2
Whole blood	26
Buffy coat	-

Whole blood (2 ml)	
--------------------	--

### Mouse / Rat

- DNA isolation for genetic test
- Genotyping

Tail, Lung, Kidney, Liver	4
Sperm	6
Tail, Lung, Kidney, Liver	29

Brain, Heart, Esophagus, Stomach, Small intestine, Large intestine, Spleen, Thymus, Lymphatic node	-
--	---

### Fish and Shellfish

- Identification of species and production region

Corbicula clam	18
Chub mackerel blood	23

Bastard halibut, Balloon fish, Ayu, Killifish, Shellfish, Loach, Eel	-
--	---

### Insects

- Genome analysis

Silkworm, Butterflies (legs), Louse	-
-------------------------------------	---

### Plants

- Identification of species and production region

Rice plant (leaf), Spinach (leaf)	32
Carnation (leaf), Peony (leaf), Camellia (leaf)	-
Cotton, Arabidopsis (leaf), Tobacco (leaf), Red bean	-
Rice kernel	-

Wheat	-
Shimeji mushroom, Pleurotaceae	-
Seaweed	-
Pufferfish (scales, muscle)	-

### Plasmid

- Amplification of target gene

E.coli	48
E.coli	49

### Fungi / Virus

- Functional analysis
- Viral DNA isolation for the identification of infector virus

SIV-infected cells	8
Branchia of KHV-infected fish	13
Yeast	20
Methicillin-resistant staphylococcus aureus (MRSA)	33
HPV-infected cells	34
Neisseria gonorrhoeae	35
Helicobacter pylori	36

Pseudomonas aeruginosa	37
Stool	38
Herpes simplex virus-type 1 (HSV-1) virus solution	39
Penicillin-resistant streptococcus pneumoniae (PRSP)	40
Vancomycin-resistant enterococcus (VRE)	41
HBV in blood serum	-

### Cell line

- Genome analysis

HepG2, Huh6 etc.	5
------------------	---

■ Numbers in the right column indicate the type of Application Guide used.

- Blue-colored numbers indicate extraction methods using the QuickGene series.
- Pink-colored numbers indicate extraction methods using the Spin Cartridge method.
- \* Where no numbers are indicated, please consult with your local contact person.

Whole blood	2
Whole blood	26
Buffy coat	-

	Title	Extraction example	Purity (A260/280)
2 ▶	Genomic DNA Isolation from Human Whole Blood	5.9 µg/200 µl Whole blood (N=5 Avg.)	1.94 (N=5 Avg.)
4 ▶	Genomic DNA Isolation from Mammalians Tissue	3.6 µg/5 mg Mouse tail (N=8 Avg.)	1.95 (N=8 Avg.)
5 ▶	Genomic DNA Isolation from Human Cultured Cell Line	5.2 µg/0.5×10 <sup>6</sup> cells HepG2	1.7
6 ▶	Genomic DNA Isolation from Mouse Sperm	3.99 µg/2.3×10 <sup>6</sup> cells	1.75
8 ▶	Viral DNA Isolation from Simian Immunodeficiency Virus (SIV) Infected Cells	7.9 µg/1×10 <sup>6</sup> cells	1.80
13 ▶	DNA Isolation from Branchia of Koi Herpes Virus (KHV) Infected Fish	1.68 µg (N=4 Avg.)	2.13 (N=4 Avg.)
16 ▶	Genomic DNA Isolation from Nail	16.7 µg/ml	1.76
17 ▶	Genomic DNA Isolation from Dental Pulp and Hard Tissues (Teeth and Bones)	37.2 µg/ml	1.87
18 ▶	DNA Isolation from Corbicula Clam	(Yield and Purity are suitable for PCR)	
20 ▶	Genomic DNA Isolation from Yeast	63.9 µg (N=7 Avg.)	1.98 (N=7 Avg.)
23 ▶	DNA Extraction from Chub Mackerel Blood Stored in TNES-6M Urea Buffer for a Long Time	12.0 µg (N=5 Avg.)	2.78 (N=5 Avg.)
26 ▶	Genomic DNA Extraction from Human Whole Blood (Spin method)	6.4 µg/200 µl Whole blood	1.97
29 ▶	Genomic DNA Extraction from Animal Tissue (Spin method)	3.5 µg/5 mg Mouse tail 4.2 µg/5 mg Mouse liver	1.93 tail 1.92 liver
32 ▶	Genomic DNA Extraction from Plants	4.3 µg/100 mg spinach (N=4 Avg.)	1.90 (N=4 Avg.)
33 ▶	Genomic DNA Extraction from Methicillin-resistant Staphylococcus aureus (MRSA)	11.6 µg/about 4 mg of wet fungi (N=3 Avg.)	1.72 (N=3 Avg.)
34 ▶	Human Papiloma Virus (HPV) DNA Extraction from Human Cervical Carcinoma Cell lines	23.5 µg HeLa	2.00 HeLa
35 ▶	Genomic DNA Extraction from Gonococcal Bacteria ( <i>Neisseria gonorrhoeae</i> )	9.0 µg/about 4.5 - 6.0mg of wet fungi (N=5 Avg.)	2.10 (N=5 Avg.)
36 ▶	Genomic DNA Extraction from <i>Helicobacter pylori</i>	3.1 µg/about 4mg of wet fungi (N=4 Avg.)	1.93 (N=4 Avg.)
37 ▶	Genomic DNA Extraction from <i>Pseudomonas aeruginosa</i>	11.4 µg	2.23
38 ▶	Bacterial Genomic DNA Extraction from Stool	16.1 µg/25 mg stool (N=2 Avg.)	2.03 (N=2 Avg.)
39 ▶	Genomic DNA Extraction from Herpes Simplex Virus-type 1 (HSV-1) Virus Solution	51 ng/about 10 <sup>7</sup> pfu/ml d41	2.14
40 ▶	Genomic DNA Extraction from Penicillin-resistant <i>Streptococcus Pneumoniae</i> (PRSP)	6.7 µg/about 4mg clinical isolate (N=2 Avg.)	1.94 (N=2 Avg.)
41 ▶	Genomic DNA Extraction from Vancomycin-resistant <i>Enterococcus</i> (VRE)	6.3 µg/about 4mg clinical isolation of wet fungi (N=2 Avg.)	1.86 (N=2 Avg.)
46 ▶	Genomic DNA Extraction from Paraffin-embedded Samples (Spin method)	1.43 µg	1.99
48 ▶	Plasmid DNA isolation from E.coli	21.4 µg	1.99
49 ▶	Plasmid DNA Extraction from E.coli (Spin method)	30.9 µg	2.00

\*Please note that some previous Application Guides have been discontinued due to outdated protocol, and the updated contents are now featured in other Application Guides.

# RNA Kit

For QuickGene-810/Mini80 [for 96 samples]

	RNA tissue kit II Processing time: 15 min/ 8 samples	RNA cultured cell kit Processing time: 17 min/ 8 samples	RNA cultured cell HC kit Processing time: 11 min/ 8 samples	RNA blood cell kit Processing time: 20 min/ 8 samples
Extraction example	ca. 100 µg/ 30 mg Mouse liver	ca. 10µg/ 1×10 <sup>6</sup> cells HL60 cell	ca. 90~150 µg/ 10 cm dish cultured HEK293 cell	ca. 4.5µg/ 1×10 <sup>7</sup> cells leukocytes

For spin-cartridge method extraction [for 96 samples]

	RNA tissue kit (spin method)	RNA cultured cell kit (spin method)	RNA cultured cell HC kit (spin method)	
Extraction example	ca. 137 µg/ 30 mg Mouse liver	ca. 10µg/ 1×10 <sup>6</sup> cells HL60 cell	ca. 213 µg/ 10 cm dish cultured HEK293 cell	

## Mammalian (Human/Cow/Poultry/Dog/Cat)

• Expression analysis such as real-time PCR and RT-PCR

Canine or feline adipose tissue, Cutis and primary-cultured adipose cells	24
Lymphatic node, Liver, Kidney	-

Leukocyte	22
-----------	----

## Mouse / Rat

• Expression analysis such as real-time PCR and RT-PCR

Liver, Brain, Lung, Kidney, Spleen, Thymus, Heart	25
Liver, Brain, Lung, Kidney, Spleen, Thymus, Heart	30
Small intestine, Esophagus, Lymphatic node, Large intestine, Stomach	-

## Insects

• Expression analysis such as real-time PCR and RT-PCR

Chironomid, Mosquito	-
----------------------	---

## Plants

• Expression analysis such as real-time PCR and RT-PCR

Wheat (leaf), Barley (leaf)	12
Arabidopsis (leaf)	-
Tomato (leaf)	-
Quinoa, Tobacco (leaf)	-
Petunia (bloom, leaf), Soybean (leaf)	-

## Cell line

• Expression analysis such as real-time PCR and RT-PCR  
• Northern blotting (cells cultured on 6 cm, 10 cm dish)  
• Microarray

Floating cell (HL60 etc.)	11
Adherent cell (COS-7, HeLa, HEK293, NIH/3T3)	14
HL60, COS-7, HeLa, HEK293, NIH/3T3	27

Cultured cell (6 cm, 10 cm dish)	21
Cultured cell (6 cm, 10 cm dish)	28

## Fungi / Virus

• The infected virus can be identified by real-time PCR and RT-PCR

VNN-infected fish	19	RS virus solution	45
Measles virus solution	42	<i>E.coli</i>	-
Influenza virus solution	43	Norovirus	-
SARS-CoV-infected cells	44		

SIV-infected cells	7
--------------------	---

■ Numbers in the right column indicate the type of Application Guide used.

- Orange-colored numbers indicate extraction methods using the QuickGene series.
- Pink-colored numbers indicate extraction methods using the Spin Cartridge method.

\* Where no numbers are indicated, please consult with your local contact person.

Liver, Brain, Lung, Kidney, Spleen, Thymus, Heart	25
Liver, Brain, Lung, Kidney, Spleen, Thymus, Heart	30
Small intestine, Esophagus, Lymphatic node, Large intestine, Stomach	-

	Title	Extraction example	Purity (A260/280)
7 ▶	Viral RNA Isolation from Simian Immunodeficiency Virus (SIV) Infected Cells	7.75 µg (N=2 Avg.)	1.84 (N=2 Avg.)
12 ▶	Total RNA Isolation from Monocotyledon and Dicotyledonous Plant Tissues	6.12 µg/50 mg Wheat leaf 12.2 µg/50 mg Barley leaf	2.11 Wheat leaf 2.12 Barley leaf
14 ▶	Total RNA Isolation from Cultured Adherent Cells (Lysing directly in culture dish)	28.1 µg/1.2×10 <sup>6</sup> cells HeLa	2.28
19 ▶	VNN (Viral Nervous Necrosis) RNA Isolation from Tilefish	(Yield and Purity are suitable for RT-PCR)	
21 ▶	Total RNA Isolation from Cultured Cells Total RNA Isolation from Cells Cultured in 6cm, 10cm Dish	129.0 µg/2.0×10 <sup>6</sup> cells HeLa	2.20
22 ▶	Total RNA Isolation from Leukocyte	6.5 µg/1.5×10 <sup>7</sup> cells Leukocyte	2.10
24 ▶	Total RNA Extraction from Canine or Feline Adipose Tissue, Cutis and Primary-cultured Adipose Cells	28.0 µg/400 mg Adipose Tissue	2.00
25 ▶	Total RNA Extraction from Various Tissues of Mouse	122 µg/30 mg Liver 21 µg/40 mg Brain	2.21 Liver 2.11 Brain
27 ▶	Total RNA Extraction from Culture Cells (Spin method)	10.5 µg/1.0×10 <sup>6</sup> cells HL60	2.11
28 ▶	Total RNA Extraction from Cultured Cells Total RNA Extraction from Cells Cultured in 6cm, 10cm Dish (Spin method)	44.0 µg/5.0×10 <sup>6</sup> cells HL60	2.17
30 ▶	Total RNA Extraction from Various Tissues of Mouse (Spin method)	137 µg/30 mg Liver	2.08
42 ▶	total RNA Extraction from Measles Virus Solution	(Yield and Purity are suitable for RT-PCR)	
43 ▶	total RNA Extraction from Influenza Virus Solution	(Yield and Purity are suitable for RT-PCR)	
44 ▶	total RNA Extraction from SARS Coronavirus (SARS-CoV) infected Cells	7.1 µg	1.90
45 ▶	total RNA Extraction from Respiratory Syncytial (RS) Virus Solution	(Yield and Purity are suitable for RT-PCR)	

\*Please note that some previous Application Guides have been discontinued due to outdated protocol, and the updated contents are now featured in other Application Guides.