

Lab Services. Your partner in liquid solutions.



This will set you free!

HERON: a unique easy plug & play robotic microplate loading system

Every laboratory that uses microplates has processes that need to be performed manually despite the installation of automated robotic systems. In particular the process of placing SBS plates onto devices requires valuable time to be committed. To address this need Lab Services has designed and developed the HERON, an instrument which features stackers, a unique plate lift and a small robotic arm with gripper that automatically introduces microplates or racks into any microplate handling device. Most importantly due to the innovative design of the HERON no stacker adjustment is needed for different rack heights. In addition the HERON is capable of handling lids, reading barcodes and disposing of microplates. With the HERON you are thus free to mix & match!

Lab Services, the specialist independent robotics solution provider with many years' experience, has developed the HERON to address a requirement of those laboratories working with microplates.

Are you ready for complete automation? Mail or call us for a liberating demonstration!
info@lab-services.nl or +31(0) 76 53 104 20

Set me free!

Why the HERON is your best manpower.

Switch operation modes easily

Easily transportable

Can be integrated with virtually any microplate device

Easy loading

Intelligent operation

Options for barcode reader, plate disposal, lid handling

Patent pending

Mail or call us for a liberating demonstration!