





Lab Services
Laboratory Automation Products




PlateButler LidPark



As specialists in automated micro plate systems and software, Lab Services has developed a custom made PlateButler™ system with the most effective software to set-up compact and complex modular robotic systems. PlateButler™ software is extremely flexible, and can be set to match your requirements, turning your lab into a 24/7 workstation.



- * Can remove all kinds of lids
- * Lids are hold by suction cups
- * Lids will be handeled with vacuum air
- * Can be integrate into Robotic Systems
- * Available as OEM product
- * Hotel can hold a maximum of 7 microplates
- * Microplates can be placed in wide or narrow



Lab Services BV, Breda, The Netherlands

Tel. +31 (0) 765310420 Fax. +31 (0) 765310421 E-mail info@lab-services.nl www.lab-services.nl