



This will turn you on!

PlateButler Turntable a revolution in plate handling.

Divers automated systems are combined from different brands of components. Lab Services encountered the many ways of plate handling within these system components. The PlateButler Turntable takes away the limitation of restricted plate orientation when designing a robotic system.

The PlateButler Turntable allows automated systems to feed plates in both landscape and portrait orientation. Thanks to the robust design the Turntable is a reliable and important part of the robotic system. Both i/o and RS232 control are possible; this allows the PlateButler Turntable to be integrated in any laboratory automation system.

Call or mail us for a dazzling demonstration!
info@lab-services.nl or +31(0)76 5310420



No limitations!

Why PlateButler Turntable keeps your laboratory running.

- Extreme reliability
- Unlimited rotation in all directions
- All SBS footprints
- Easy robot handling
- OEM availability
- Robust metal construction
- Small footprint
- RS232 control and i/o control