


PRACTICAL INNOVATION

Ziath

# DATA PAQ™ HIGH SPEED SINGLE RACK SCANNER

Product ID – ZTS-A6



The A6-sized scanner is currently our most popular product. It offers easy set up, straight out of the box, and with only a few minutes needed for configuration, you'll be ready to start scanning in no time at all. A one-time automatic calibration means that you'll then be set up and ready to read any and all tubes – even those hard to read REMP tubes.

Our A6 is fully integratable into robotics platforms.  
Base plates for platforms such as the TECAN Evo are available.

[www.ziath.com](http://www.ziath.com)

## System Specification


---

### High Speed Single Rack Scanner (A6)

- Operating Systems: Windows XP/Vista
- Scanning Device: Modified Fujitsu flatbed scanner
- Compatible rack vendors
  - Thermo Matrix
  - Thermo Abgene
  - Thermo Nunc
  - Remp
  - Micronic
  - FluidX
  - .... And all other 2D coded tubes
- Compatible rack sizes: 384, 96, 48, 24
- Export data to Excel, text, XML, rack images
- Fast Scanning — 5 seconds per rack
- Displays damaged/missing tubes
- One-click scan for any size or brand of tube
- Simple installer — no dongles required
- Multiple scanners/multiple configurations
- 110/220 v automatically selected
- Fully integratable into robotics platforms

## System Footprint

---



**For Base Plate Dimensions for Specific Integration Applications, Please Inquire.**

**For more information please contact Ziath Ltd.**

[www.ziath.com](http://www.ziath.com) • [sales@ziath.com](mailto:sales@ziath.com)

**T +44 (0) 1223 968 125 • F +44 (0) 1223 771 615**