


PRACTICAL INNOVATION



EFLY2 MANUAL HEAT SEALER

Product ID – ZEHS-01




Exclusive!
Only
available
from
Ziath!

www.ziath.com

System Specification

- Desktop semi-automatic Heat Applied Plate Sealing System
- Seal SBS format micro well plates
- Seal plate height with no adjustment from 9 – 48mm high
- Temperature range from ambient to 200° C.
- Seal time from 0.1 – 9.9 seconds
- Simple slide operation
- Automatic eject when timed seal completed
- Combined temperature and seal time display
- Light weight system only 6kgs (14lbs)
- 110/220 vac power compatible
- No compressed air required

System Footprint



For more information please contact Ziath Ltd.

www.ziath.com • sales@ziath.com

Tel +44 (0) 1223 968 125 • Fax +44 (0) 1223 771 615