


DATAPAQ™ HIGH SPEED 1D RACK BARCODE SCANNER

Product ID – ZTS-1DR



The Ziath 1D Rack Barcode Scanner is designed to be used in conjunction with our popular A6-sized scanner. The Rack Barcode Scanner can be attached to any side of the A6 scanner and is fully adjustable in terms of the read height and angle to accommodate different placement of barcodes on racks. Alternatively, the scanner can be attached permanently to the A6.

The Rack Barcode Scanner offers easy set up, straightforward USB connectivity, and seamless integration with the DataPaq software employed by the A6 scanner.


System Specification

High Speed 1D Rack Barcode Scanner

- Operating Systems: Windows XP/Vista
- Scanning Device: Opticon Rack Barcode Scanner
- Use in conjunction with the Ziath A6 High Speed Single Rack Scanner
- Fully integrated with DataPaq software
- Via DataPaq, rack and tube data can be exported to Excel, text, XML, rack images
- Fast Scanning — <1 second per rack barcode
- Simple installation via USB

System Footprint

Note: Dimensions shown are approximate. Actual footprint/height may vary due location and angle required to read specific rack barcodes.



For more information please contact Ziath Ltd.

www.ziath.com • sales@ziath.com

Tel +44 (0) 1223 968 125 • Fax +44 (0) 1223 771 615