

# DATA PAQ™

## HIGH SPEED SINGLE TUBE SCANNER

Product ID – ZTS-S1



Our Single-Tube Scanner, like all Ziath scanners, offers easy set up, straight out of the box – along with just the kind of portability some lab work requires. With only a few minutes needed for installation, you'll be ready to start scanning in no time at all.

Ideal for benchtop use, or simply take it with you along with your laptop. Like all our scanner products, the Single Tube scanner is very fast – just hold the tube in the lens window and the read shows immediately on screen.

## System Specification


---

### DataPaq High Speed Single Tube Scanner (S1)

- Operating Systems: Windows XP/Vista
- Compatible rack vendors
  - Thermo Matrix
  - Thermo Abgene
  - Thermo Nunc
  - Micronic
  - FluidX
  - Matrical
- Export data to Excel, text, tube image
- Fast Scanning — <1 second per tube
- Null read for damaged tube
- Instant read for any size or brand of tube
- Simple installer — no dongles required
- USB connectivity – no independent power supply needed.

## System Footprint

---



For more information please contact Ziath Ltd.

[www.ziath.com](http://www.ziath.com) • [sales@ziath.com](mailto:sales@ziath.com)

T +44 (0) 1223 968 125 • F +44 (0) 1223 771 615