



PiresTrack™

The PiresTrack™ line is unique for being the sole line of affordable pipette reservoirs that can come with lids.

It consists of both sterile and non-sterile pipette reservoirs, that are available in 25 - 50 ml and 50 - 100 ml sizes and 8- or 12channel liner reservoirs.

Non-sterile bulk packs are also available. For the products without lids, lids can be ordered separately.

The lids protect your reagent from foreign matter falling in, and offer protection against spilling when filled reservoirs need to be moved. They also make it possible to stack the reservoirs (like in the picture above).

All our open reservoirs have ridges on one side that make it easier to take off drops hanging on the pipette tips. In addition, they all have a non-nesting design, which means that the single reservoirs can easily be separated from the bulk pile.

The PiresTrack™ products are made of high impact polystyrene, a material that has proven it's suitability for this kind of product.

On request other type material can be made available. Summarizing, PiresTrack™ offers the following advantages :

- unique lid design which makes reservoirs stackable
- prevents contamination and spilling
- ridges for convenient accurate pipetting
- non-nesting design
- individually packed sterile
- bulk packed non-sterile

PiresTrack™ Pipette reservoirs for single and multi-channel pipettes


Catalog No.	6110	6111	6010
Model	open	open	open
With Lid	no	yes	no
Reservoir volume	50ml	50ml	50ml
Stackable	no	yes	no
Material	HIPS	HIPS	HIPS
Sterile	yes	yes	no
Indiv./pkg.	yes	yes	10/pkg.
Units/cs.	100	100	200



PiresTrack™ Pipette Reservoirs

Catalog No.	6112	6113	6011	6114	6115	6116	6117
Model	open	open	open	8 channel	12 channel	8 channel	12 channel
With Lid	no	yes	no	no	no	yes	yes
Reservoir volume	100ml	100ml	100ml	5ml/ch.	4ml/ch.	5ml/ch.	4ml/ch.
Stackable	no	yes	no	no	no	yes	yes
Material	HIPS	HIPS	HIPS	HIPS	HIPS	HIPS	HIPS
Sterile	yes	yes	no	yes	yes	yes	yes
Indiv./pkg.	yes	yes	10/pkg.	yes	yes	yes	yes
Units/cs.	100	100	200	100	100	100	100

Catalog No.	6014	6015	6030	6031
Model	8 channel	12 channel	reservoir lid	reservoir lid
With Lid	no	no	na	na
Reservoir volume	5ml/ch.	4ml/ch.	25-50ml	50-100ml / 8-12ch
Stackable	no	no	na	na
Material	HIPS	HIPS	HIPS	HIPS
Sterile	bulk	bulk	bulk	bulk
Indiv./pkg.	10	10	10	10
Units/cs.	100	100	50	50



For more in depth information on these products please consult our Internet pages on www.concenttrack.com or contact us to provide you with our product brochure. info@concenttrack.nl / Tel. +31(0)365329011

HIPS = high impact polystyrene